



Easy kit for 100 mL

According to US-EPA 1642 and 1643 method.
Raw and treated wastewater, surface water, drinking water, recreational water, shellfish extracts, sediments, and sludge

Sample Volume 100 mL 10 assays
2x TSA ready to use
Host strain and positive control included

BP1609 Easy kit for somatic coliphages

BP1619 Easy kit for F-specific coliphages



Biological material

BP1629:
Host strain *E. coli* CN13 20 units
For analysis of Somatic coliphages

BP1631
Host strain *E. coli* HS 20 units
For analysis of F-specific coliphages

BP1626:
Virus bacteriophage ϕ X174 20 units
Somatic coliphage

BP1627:
Virus bacteriophage MS2 20 units
F-specific coliphage

BP1608:
Host strain *E. coli* CN13 10 units
Virus bacteriophage ϕ X174 10 units
For analysis of Somatic coliphages

BP1618:
Host strain *E. coli* CN13 10 units
Virus bacteriophage MS2 10 units
For analysis of F-specific coliphages

