

Labcon

Redefining Laboratory Performance

LOW TEMPERATURE INCUBATORS

High-performance engineering
for results you can trust



LABORATORY INCUBATORS

FEATURES

- 304 stainless steel inner chamber with powder coated mild steel outer cabinet and non-tip shelves
- Inner glass door with full peripheral silicon rubber gasket
- Outer glass door with peripheral silicon rubber gasket
- Hermetically sealed compressor with overpad protection
- Microprocessor temperature control with digital temperature indicator
- Over / under temperature alarms
- Auto / manual defrost system
- Safety thermostat (DIN 12880 standard)

Note: This range of incubators should be used in preference to the FSIM range of incubators when the equipment is required to operate below ambient temperature.

Applications: Biology; Botany; Agriculture; Food Industry, etc

Control Accuracy: $\pm 0.2^{\circ}\text{C}$ @ 37°C

Uniformity: $\pm 0.5^{\circ}\text{C}$ @ 37°C

OPTIONS


1. Illumination
2. Special requirements on request



THREADbridge
www.kreon.co.za/threadbridge

REAL-TIME DATA INTELLIGENCE

Q1 YEAR WARRANTY*
QUALITY GUARANTEE



Pictures may vary from actual product.


Technical Data

MODEL	DIMENSIONS						VOLUME (APPROX. ℓ)	SHELVES		TEMPERATURE RANGE °C	POWER		WEIGHT
	EXTERNAL			INTERAL				MAX	SUPPLIED		VOLTAGE	Kw	
	H	W	D	H	W	D							
LTIM 10	940	615	660	510	480	410	100ℓ	5	2	0 ~ +80°C	220/240	0,5	86
LTIM 20	1275	750	760	740	625	500	230ℓ	6	2	-5 ~ +80°C	220/240	1	140
LTIM 40	1775	750	760	1240	625	500	400ℓ	12	2	-5 ~ +80°C	220/240	1,5	180
LTIM 70*	1775	1515	760	1240	1400	500	850ℓ	12	2	0 ~ +80°C	220/240	2	290

LTIM 70* - Double door unit, no center door column

DISCLAIMER

Output performance figures quoted are dependent on various factors. Scientific reserves the right, without notice, to make changes and revisions to product specifications, materials and design, when we believe it will provide better performance, durability, and efficiency. Unit measurements may vary depending on point of reference.



*Labcon Equipment comes with a standard ONE year warranty.

Wireless
Sensor
Network,
for
Wired
Devices

